












Medical Device Symbols included on Serim Product labels and Instructions for Use (IFU)

SYMBOL	SYMBOL TITLE	EXPLANATORY TEXT	STANDARD TITLE	STANDARD REFERENCE
	Manufacturer	Indicates the medical device manufacturer.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.1 FDA Recognition # 5-117 ISO 7000 Reference #3082 FDA Recognition # 5-103
	Use-by Date	Indicates the date after which the medical device is not to be used. Date format is YYYY-MM-XX	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.4 FDA Recognition # 5-117 ISO 7000 Reference #2607 FDA Recognition # 5-103
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.,	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.5 FDA Recognition # 5-117 ISO 7000 Reference #2492 FDA Recognition # 5-103
	Catalog #	Indicates the manufacturer's catalog # so that the medical device can be identified.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.1.6 FDA Recognition # 5-117 ISO 7000 Reference #2493 FDA Recognition # 5-103
	Keep away from sunlight	Indicates a medical device that needs protection from light sources.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.2 FDA Recognition # 5-117 ISO 7000 Reference #0624 FDA Recognition # 5-103
	Keep dry	Indicates a medical device that needs to be protected from moisture.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.4 FDA Recognition # 5-117 ISO 7000 Reference #0626 FDA Recognition # 5-103
	Temperature limit	Indicates the temperature limits to which the medical device can be safely exposed.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.3.7 FDA Recognition # 5-117 ISO 7000 Reference #0632 FDA Recognition # 5-103
	Do not re-use	Indicates the temperature limits to which the medical device can be safely exposed.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.4.2 FDA Recognition # 5-117 ISO 7000 Reference #1051 FDA Recognition # 5-103
	Consult instructions for use	Indicates the need for the user to consult the instructions for use.	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.4.3 FDA Recognition # 5-117 ISO 7000 Reference #1641 FDA Recognition # 5-103
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot,	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements	ISO 15223-1 Reference #5.4.4 FDA Recognition # 5-117 ISO 7000 Reference #0434A FDA Recognition # 5-103
	Medical device	Indicates that the item is a medical device	Medical Devices — Symbols to be used with medical device labels, labeling, and information to be supplied — Part 1: General requirements. (Draft of New Version)	ISO/DIS 15223- 1:2020(E) DRAFT Reference no. 5.7.7